

# Introduction

---

This document outlines a step-by-step workflow to process Collecta Clonetracker Barcode single-cell data and assign clonetracker barcodes to each cell using shell and Python scripts. The aim is to ensure reproducibility and clarity in sequencing data processing and barcode analysis.

## Input Data

---

The input data for this workflow includes:

- **FASTQ files from Gene Expression Profiling NGS data:**
  - /mnt/project/Collecta\_scCRISPR/VCU\_scCloneTracker/data/AD\_GEX/AD\_GEX\_S3\_R1\_001.fastq.gz
  - /mnt/project/Collecta\_scCRISPR/VCU\_scCloneTracker/data/AD\_GEX/AD\_GEX\_S3\_R2\_001.fastq.gz
- **FASTQ files from the Collecta Clonetracker barcode NGS library:**
  - /mnt/project/Collecta\_scCRISPR/VCU\_scCloneTracker/data/AD\_FBP1\_S7\_R1\_001.fastq.gz
  - /mnt/project/Collecta\_scCRISPR/VCU\_scCloneTracker/data/AD\_FBP1\_S7\_R2\_001.fastq.gz
- **Reference barcode sequences provided by Collecta:**
  - Collecta-CloneTrackerXP-5M-Pool1-BC14-LNGS-300-Library-Design.txt
  - Collecta-CloneTrackerXP-50M-BC30-LNGS-300-Library-Design.txt

## Step 1: Load Reference Barcode Sequences

---

Analyze Gene Expression Profiling NGS data and extract cell barcodes using `cellranger count`:

```
/mnt/project/Pipeline/Software/cellranger-8.0.0/cellranger count \  
--id=AD_GEX \  
--fastqs=/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/data/AD_GEX/ \  
--transcriptome=/mnt/project/Pipeline/Reference/10XGenomics/refdata-gex-  
GRCh38-2024-A \  
--create-bam true \  
--include-introns false  
  
zcat  
/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/AD_GEX/outs/filtered_feature_bc_  
matrix/barcodes.tsv.gz > AD_barcode.xls  
  
sed 's/-.*/' AD_barcode.xls > AD_barcode_cleaned.tsv
```

## Step 2: Process Barcode NGS Library Using UMI-Tools

---

Use `umi_tools` to extract and match barcodes from FASTQ files:

## Shell Script

```
#!/bin/bash

FILEIN1=/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/data/AD_FBP1_S7_R1_001.fastq.gz
FILEOUT1=/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/data/`basename ${FILEIN1}`_extracted.fastq.gz
FILEIN2=/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/data/AD_FBP1_S7_R2_001.fastq.gz
FILEOUT2=/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/data/`basename ${FILEIN2}`_extracted.fastq.gz

WHITELIST=/mnt/project/Collecta_scCRISPR/VCU_scCloneTracker/data/AD_barcode_cleaned.tsv

umi_tools extract \
  --bc-pattern=CCCCCCCCCCCCCCCCNNNNNNNNNNNN \
  --stdin $FILEIN1 \
  --stdout $FILEOUT1 \
  --read2-in $FILEIN2 \
  --read2-out=$FILEOUT2 \
  --whitelist=$WHITELIST
```

## Step 3: Identify the Best Sequence for Each UMI

---

Run the following Python script to select the best-quality sequence for each UMI:

```
python best_sequence_umi.py
```

## Step 4: Assign Clonetracker Barcodes to 10X Cell Barcodes

---

Use this script to map Clonetracker barcodes to cell barcodes:

```
python barcode_process_umi_5.py
```

## Step 5: Final Barcode Assignment Using UMI Distribution

---

Perform the final barcode assignment based on UMI distribution and criteria:

```
python umi_distribution_analysis_umi_5.py
```

## Conclusion

---

This workflow provides a systematic approach to process Collecta Clonetracker single-cell data, assign barcodes, and analyze barcode distributions. Scripts and parameters can be adjusted based on dataset-specific requirements. By following these steps, you can ensure accurate and reproducible results for your analysis.